

# AL WASAY

## Diagnostic Services

---

## Welcome to Al Wasay Diagnostic Services

Office No. 4, Street No. 02, Opp Wapda  
Town Phase II, Northern Bypass, Pace City,  
Multan.






0333-7349146 – Muhammad Jahangir Bhatti



0307-5485849 – Muhammad Kamran Arshad

## Blood Roller Mixer – 5 Tubes Capacity

-  Reliable Laboratory Solutions
-  For Medical Laboratories & Blood Banks
-  Gentle & Continuous Blood Mixing







**Al Wasay Diagnostic Services**  
Office No. 4, Street No. 02, Opp Wapda Town Phase II,  
Northern Bypass, Pace City, Multan

 WhatsApp: 0333-7349146 - Muhammad Jahangir Bhatti  
 Call: 0307-5485849 - Muhammad Kamran Arshad

# Product Overview

Designed for safe, reliable, and continuous blood sample handling.

-  Designed for gentle and continuous blood mixing.
-  Prevents clot formation during handling.
-  Ensures uniform and consistent sample preparation.
-  Ideal for routine laboratory and blood bank use.



### Al Wasay Diagnostic Services

Office No. 4, Street No. 02, Opp Wapda Town Phase II,  
Northern Bypass, Pace City, Multan

 **WhatsApp:** 0333-7349146 – Muhammad Jahangir Bhatti

 **Call:** 0307-5485849 – Muhammad Kamran Arshad

## Why Blood Roller Mixer?

Key benefits for safe, accurate, and reliable blood handling



Prevents blood **clot formation** during sample handling



Ensures **consistent** and **uniform mixing** results



**Protects** blood sample **integrity** for accurate testing



**Reduces** manual handling and human error



**Smooth, gentle rotation** safe for sensitive samples









Ideal for **continuous** routine laboratory and blood bank operations



# U-Tech Biomedical

Reliable Laboratory Solutions

## Technical Specifications.

-  **Tube Capacity:** 5 Tubes
-  **Mixing Method:** Continuous Roller Mixing
-  **Roller Length:** 115 mm
-  **Roller Material:** High-quality, chemical-resistant rubber rollers
-  **Speed Control:** Fixed optimized speed for blood sample safety
-  **Operation Type:** Continuous & smooth rotation



### Al Wasay Diagnostic Services

Office No. 4, Street No. 02, Opp Wapda Town Phase II,  
Northern Bypass, Pace City, Multan

 **WhatsApp:** 0333-7349146 – Muhammad Jahangir Bhatti

 **Call:** 0307-5485849 – Muhammad Kamran Arshad

## Electrical Specifications



Power Supply: 220V AC, 50 Hz



Power Consumption: 12V DC



Current: 2.5 Amp



Energy Efficient Design



### Al Wasay Diagnostic Services

Office No. 4, Street No. 02, Opp Wapda Town Phase II,  
Northern Bypass, Pace City, Multan

**WhatsApp:** 0333-7349146 – Muhammad Jahangir Bhatti

**Call:** 0307-5485849 – Muhammad Kamran Arshad

## Physical Specifications



Length:  
**333 mm**



Width:  
**323 mm**



Height:  
**115 mm**



Net Weight:  
**2.775 kg**



Gross Weight:  
**3.960 kg**



Design:  
**Compact &  
space-saving**



## Build & Safety Features



Durable ABS body for long service life



Low vibration & silent operation for lab-friendly use



Heat-resistant motor suitable for continuous operation



Stable base design for safe and secure placement



Maintenance-free performance reduces downtime



## Applications



**Clinical Laboratories**  
routine analysis &  
testing



**Hospital Laboratories**  
in-patient & out-patient  
testing



**Blood Banks**  
sample mixing &  
preservation







**Diagnostic Centers**  
specialized diagnostic  
procedures








## Optional Features

### Battery Backup Model (Optional)

-  Ensures uninterrupted operation during power failures
-  Ideal for areas with frequent power fluctuations
-  Maintains continuous and safe blood mixing
-  Highly recommended for blood banks and critical labs



## After-Sales Support & Service

-  ISO-certified manufacturing support
-  1-Year Warranty for complete peace of mind
-  Fast after-sales service within 24 hours
-  Dedicated technical support team
-  Reliable spare parts & long-term service availability



# Thank You

## For Choosing Reliable Laboratory Solutions

Committed to Accuracy, Safety & Trust in Healthcare