




Authorized Supplier & Diagnostic Equipment Provider

A soft tagline is seart professional medical diagnostic company.

AL WASAY DIAGNOSTIC SERVICES

Office No.4 Street No.02 Opp Wapda Town Phase
II Northern Bypass Pace City Multan

Muhammad Jahangir Bhatti – 03337349146 

Muhammad Kamran Arshad – 03075485849 



AL WASAY DIAGNOSTIC SERVICES

Hemasure HS-3000

Accuracy you can trust.

High-performance automated hematology analyzer designed for accuracy, efficiency, and routine laboratory use.



AL WASAY DIAGNOSTIC SERVICES

Office No.4 Street No.02 Opp Wapda Town Phase II, Northern Bypass Pace City Multan

Muhammad Jahangir Bhatti – 03337349146



Muhammad Kamran Arshad – 03075485849

Hemasure HS-3000

Accuracy you can trust.

User-friendly operation



Compact design



Large data storage capacity



Flexible aspiration mode



No Sleep Mode – Continuous operation without standby delay



Editable Results – User can edit and correct test results



Low Reagent Consumption – Saves cost with minimal reagent use



High Throughput Efficiency – 1 reagent pack supports 1000+ cycles

AL WASAY DIAGNOSTIC SERVICES

Office No.4 Street No.02 Opp Wapda Town Phase II, Northern Bypass Pace City Multan

Muhammad Jahangir Bhatti – 03337349146  Muhammad Kamran Arshad – 03075485849 

Hemasure HS-3000

Accuracy you can trust.



Reagents: 2



Reportable Parameters: 21



Throughput: 60 Samples Per Hour

Optimized for routine laboratory workflow with reliable and consistent performance.

AL WASAY DIAGNOSTIC SERVICES

Office No.4 Street No.02 Opp Wapda Town Phase II Northern Bypass Pace City Multan
Muhammad Jahangir Bhatti – 03337349146  Muhammad Kamran Arshad – 03075485849 

Hemasure HS-3000

Accuracy you can trust.



Applications



Clinical
Laboratories



Hospitals



Diagnostic
Centers



Medical Testing
Facilities

Designed to meet the daily diagnostic needs of modern healthcare facilities.



AL WASAY DIAGNOSTIC SERVICES

Hemasure HS-3000

Accuracy you can trust.



Performance & Reliability

✓  High Accuracy & Consistent Results

✓  Reliable Performance for Routine Laboratory Use

✓  Designed for Continuous Operation

✓  Cost-Efficient with Low Reagent Consumption

✓  Engineered for Modern Diagnostic Workflows

Built to deliver dependable results with confidence in every test.

AL WASAY DIAGNOSTIC SERVICES

Office No.4 Street No.02 Opp Wapda Town Phase II Northern Bypass Pace City Multan

Muhammad Jahangir Bhatti – 03337349146 

Muhammad Kamran Arshad – 03075485849 

Trusted Diagnostic Solutions

Powered by Innovation. Driven by Accuracy.

Hemusure HS-3000

Accuracy you can trust.



- ✓ Proven diagnostic performance
- ✓ Trusted by laboratories and healthcare professionals
- ✓ Designed for reliability and long-term use



“Contact us today to learn more about diagnostic solutions.”




THANK YOU FOR VISITING US

We look forward to serving your diagnostic needs with trust and quality.

AL WASAY DIAGNOSTIC SERVICES

Office No.4 Street No.02 Opp Wapda Town Phase II
Northern Bypass Pace City Multan

Muhammad Jahangir Bhatti – 03337349146 

Muhammad Kamran Arshad – 03075485849 

Vector rectangles in brand colors