

Al Wasay Diagnostic Services

Authorized Supplier & Diagnostic Equipment Provider

DIGITAL WATER BATH

U-TECH Bio-Medical



AWD
AL WASAY
DIAGNOSTIC SERVICES



0333-7349146
محمد جهانگیر بھٹی



0307-5485849
محمد کامران ارشد

Al Wasay Diagnostic Services

Office No. 4, Street No. 02, Opp Wapda Town Phase II, Northern Bypass, Pace City, Multan

Al Wasay Diagnostic Services

Authorized Supplier & Diagnostic Equipment Provider

Product Overview

Digital Water Bath – U-TECH Bio-Medical

What is a Digital Water Bath?

A Digital Water Bath is a laboratory heating device designed to maintain samples at a constant and controlled temperature. It uses water as a heat transfer medium to provide uniform and accurate heating. Digital water baths are widely used in medical, diagnostic, and research laboratories for incubation, sample warming, reagent preparation, and temperature-sensitive testing procedures. The digital temperature control ensures precision, safety, and repeatable results in laboratory workflows.



Al Wasay Diagnostic Services

Authorized Supplier & Diagnostic Equipment Provider

Key Features

Digital Water Bath – U-TECH Bio-Medical



Uniform & Accurate Heating

Provides consistent temperature distribution for reliable laboratory results.



Digital Temperature Control

Ensures precise monitoring and accurate temperature regulation.



Efficient Water Circulation System

Maintains uniform heat transfer across the entire bath chamber.



Energy Efficient Operation

Optimized heating system reduces power consumption.



Corrosion-Resistant Stainless Steel Body

Designed for long-term use in laboratory environments.



Compact & Durable Design

Space-saving construction with high operational reliability.



Al Wasay Diagnostic Services

Authorized Supplier & Diagnostic Equipment Provider

Technical Specifications

Digital Water Bath – U-TECH Bio-Medical

Parameter	Specification
Temperature Range	Ambient to 100°C
Temperature Accuracy	±0.1°C
Temperature Control	Digital PID Controller
Display Type	LED / Digital Display
Bath Capacity	As per model
Heating System	High-efficiency immersion heater
Body Material	Corrosion-resistant stainless steel
Power Supply	220-240V / 50-60Hz
Safety Features	Over-temperature protection
Application	Clinical, Diagnostic & Research Laboratories



Applications

Digital Water Bath – U-TECH Bio-Medical



Clinical Laboratories

Used for controlled heating of samples, reagents, and test solutions in routine lab procedures.



Diagnostic Centers

Ideal for incubation and temperature-sensitive diagnostic testing processes.



Blood Banks

Suitable for warming blood samples and maintaining required temperatures.



Research & Educational Laboratories

Used in research experiments, training labs, and academic institutions.



Pharmaceutical Laboratories

Supports preparation and testing of pharmaceutical samples under controlled conditions.



Hospitals & Medical Institutions

Used in pathology, microbiology, and general laboratory departments.



0333-7349146

محمد جہانگیر بھٹی



0307-5485849

محمد کامران ارشد



Company Name:
AI Wasay Diagnostic Services



Address:
Office No. 4, Street No. 02,
Opp Wapda Town Phase II,
Northern Bypass, Pace City, Multan

Al Wasay Diagnostic Services

Authorized Supplier & Diagnostic Equipment Provider

Why Choose Us

Al Wasay Diagnostic Services



Trusted Diagnostic Equipment Supplier

Recognized for supplying reliable and high-quality medical and laboratory equipment.



Quality Assurance & Certified Products

All products meet international laboratory and medical standards.



After-Sales Support & Technical Assistance

Dedicated support team for installation, guidance, and maintenance.



Warranty & Spare Parts Availability

Ensures long-term peace of mind and uninterrupted laboratory operations.




Competitive Pricing

Cost-effective solutions without compromising on quality.



Customer-Focused Approach

Committed to long-term partnerships with hospitals and diagnostic centers.



0333-7349146
محمد جهانگیر بھٹی



0307-5485849
محمد کامران ارشد

Al Wasay Diagnostic Services

Office No. 4, Street No. 02, Opp Wapda Town Phase II,
Northern Bypass, Pace City, Multan

Al Wasay Diagnostic Services

Authorized Supplier & Diagnostic Equipment Provider

Thank You

For your time and interest in our diagnostic solutions
We look forward to serving your laboratory and diagnostic needs

 **0333-7349146** محمد جہانگیر بھٹی

 **0307-5485849** محمد کامران ارشد



Al Wasay Diagnostic Services

Office No. 4, Street No. 02, Opp Wapda Town Phase II,
Northern Bypass, Pace City, Multan