




XD690 **Electrolyte** **Analyzer Series**



Precision. Stability. Reliability.

About XD690 Electrolyte Analyzer Series

The XD690 Series is a high-precision electrolyte analyzer designed for reliable and stable laboratory diagnostics. It ensures accurate measurement of K^+ , Na^+ , Cl^- , Li^+ , iCa^{++} , pH and TCO_2 depending on model configuration. The system provides fast testing speed, stable results and user-friendly interface for modern laboratories.



MAIN FEATURES



Computer monitored testing
for high precision & stability



Automatic air bubble
detection



Large HD touch screen
display



Eco-friendly sealed
reagent pack



Automatic 24-hour
operation



Automatic data saving
during power failure



Self-diagnosis & electrode
alarm system



QC adjustment with
31-day QC charts

Technical Specifications

Specification	Details
Sample Type	Plasma, Serum, Blood
Test Speed	≥ 60 samples/hour
Sample Volume	≥ 120 µl
Display	Large LCD Touch Screen
Output	Built-in
Power	RS232C
Environment	AC 220±22V, 50Hz
Size	5–40°C, Humidity ≤80%
Weight	390 × 270 × 408 mm

Applications

Suitable for:

- Hospitals
- Diagnostic Laboratories
- Emergency Units
- Research Centers
- Clinical Laboratories



Available Models & Configurations

XD690 Series – Multi-Parameter Electrolyte Analysis Solutions

• XD690 I

K^+ / Na^+ / Cl^- Analyzer

• XD690 II

K^+ / Na^+ / Cl^- / Li^+ Analyzer

• XD690 III


K^+ / Na^+ / Cl^- / iCa^{++} / pH
Analyzer

• XD690 IV

K^+ / Na^+ / Cl^- / iCa^{++} /
 Li^+ / pH Analyzer

• XD697

K^+ / Na^+ / Cl^- / TCO_2
Analyzer



Contact Us

AWD – Al Wasay Diagnostic Services

“Committed to Precision Diagnostics”

 Address:

Office No. 4, Street No. 02, Opp Wapda Town Phase II,
Northern Bypass, Pace City, Multan

 Sales & Support Contacts:

0333-7349146 – Muhammad Jahangir Bhatti 

0307-5485849 – Muhammad Kamran Arshad 